

# 4-gel vertical electrophoresis packages

Whether you require a package to match the market leader or one to address any budget or vertical electrophoresis technique, Cleaver Scientific now offers the ultimate range of Combination Packages that will meet your needs.

## Run 4 gels in the tank while casting another 4 externally



### Ordering Information

#### CVS10TETRAD1 OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels, includes: -

- CVS10DSYS-CU – 10x10cm mini-vertical unit with casting base and external casting stand, plus 2x 1mm 12-sample combs, 2x plain glass plates with 1mm spacer and 2x notched plates (Page 22)
- VS10-12-1 – 2 additional 1mm 12-sample combs
- VS10PGS1 – 1 pack of 2x plain glass plates with 1mm spacer
- VS10NG – 1 pack of 2x notched glass plates
- VS10NGS1 – 2 packs of 2x notched glass plates with 1mm spacer

For 0.75, 1.5 and 2mm versions use CVS10TETRAD0.75, CVS10TETRAD1.5 and CVS10TETRAD2 respectively

## With interchangeable 4-blot module




### Ordering Information

#### CVS10TETRAD1CBS OmniPAGE TETRAD with interchangeable 4-blot module, includes: -

- CVS10TETRAD1 – OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels
- VS10BI – 4-blot insert with 4 blotting cassettes and 8 fibre pads (Page 34)

For 0.75, 1.5 and 2mm versions use CVS10TETRAD0.75CBS, CVS10TETRAD1.5CBS and CVS10TETRAD2CBS respectively

## With standard Midi power supply option



### Ordering Information

#### CVS10TETRAD1-CS300 OmniPAGE TETRAD with standard Midi power supply option, includes: -

- CVS10TETRAD1 – OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels
- CS-300V – 300V, 700mA, 150W midi power supply (Page 39)

For 0.75, 1.5 and 2mm versions use CVS10TETRAD0.75-CS300, CVS10TETRAD1.5-CS300 and CVS10TETRAD2-CS300 respectively

## With standard Midi power supply option and interchangeable 4-blot module



### Ordering Information

#### CVS10TETRAD1CBS-CS300 OmniPAGE TETRAD with Midi power supply and interchangeable 4-blot module, includes: -

- CVS10TETRAD1CBS – OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels with blotting module
- CS-300V – 300V, 700mA, 150W midi power supply (Page 39)

For 0.75, 1.5 and 2mm versions use CVS10TETRAD0.75CBS-CS300, CVS10TETRAD1.5CBS-CS300 and CVS10TETRAD2CBS-CS300 respectively

## With high current Maxi power supply option



### Ordering Information

#### CVS10TETRAD1-CS3AMP OmniPAGE TETRAD with high current Maxi power supply option, includes: -

- CVS10TETRAD1 – OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels
- CS-3AMP – 300V, 3000mA, 300W maxi power supply (Page 39)

For 0.75, 1.5 and 2mm versions use CVS10TETRAD0.75-CS3AMP, CVS10TETRAD1.5-CS3AMP and CVS10TETRAD2-CS3AMP respectively

## With high current Maxi power supply option and interchangeable 4-blot module



### Ordering Information

#### CVS10TETRAD1CBS-CS3AMP OmniPAGE TETRAD with Maxi power supply and interchangeable 4-blot module, includes: -

- CVS10TETRAD1CBS – OmniPAGE TETRAD Mini-Protein Electrophoresis Package for 1mm self-cast gels with blotting module
- CS-3AMP – 300V, 3000mA, 300W maxi power supply (Page 39)

For 0.75, 1.5 and 2mm versions please to remember to use the correct number. E.g. - CVS10TETRAD0.75CBS-CS3AMP for a 0.75mm version